
# PRIXTON

POWER FLEX

**PISTOLA DE MASAJE** MASSAGE GUN

---

**INDEX**

ES	3
<hr/>	
EN	13
<hr/>	
FR	24
<hr/>	
IT	33
<hr/>	
DE	42
<hr/>	
PT	51
<hr/>	
NL	60
<hr/>	

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
[www.prixton.com](http://www.prixton.com) LA TRASTIENDA DIGITAL

**WARNING To charge this device, make sure to use a power adapter with the following values: 5 volts and a recommended amperage of at least 1000mA.**

## Instructions

### Power operation cautions

Press the power switch at the top of the machine and hold 1 second to turn it on and start working.(Power light on) Press the power switch at the top of the machine and hold 1 second again to turn it off.(Power light off)Press the heat switch on the top of the machine to turn it on and start the device (key indicator light on). Press the heat switch at the top of the machine again to turn off the device (key indicator off).

### Device

Equipment Power switch on,quick press the power bottom once to adjust it to level 1,press twice to adjust it to level 2, and so on,a total of 6 adjustable,display 6 levels.Level 1-6 corresponding frequency levels.Heating button,press once to level 1,and press level 2,and level 3.When the heat working,the temperature of the three gears are 40 degrees, 50 degrees, 60 degrees, the heating time to 40 degrees is about 2 minutes.

### Charging

Before the first use,charging 4h to fill it up,charging the device does not work, the power is divided into 4 levels. The battery can be charged at any time and at any charge level.It is not recommended that the battery be completely depleted.Make sure the power is off when the machine is not in use.The operating time of the device depends on the speed level and the pressure applied during use.


U-shaped head for spinal area

Small ball type head for small muscle groups

Small round head for joint area


Small flat head for all parts of the body

Large ball-shaped head for large muscle groups


### Use scenes


1. Back massage

2. Neck massage



3. Foot massage

4. Arm massage

### Adjustable pulling belt


Pull the strap to adjust the installation adjust the installation.

Curved head  
for bending parts

Strap head  
for fixing pulling belt

Massage straps  
for other parts of the neck,  
back, etc

Massage block  
for all parts of the body

### Configuration and Specifications

#### Configuration

1. Control panel
2. Heated electrode piece
3. Removable and replaceable vibration head. When disassembling the vibration head, pull it out hard
4. DC charging port

#### Product parameters

Battery voltage: 7.4V  
Rate power: 24W  
Rate voltage: 5V  
Battery type: Lithium battery

## Physical Therapy Precautions

### Post-injury

The first physical therapy session after injury should be gradual according to the patient's satisfaction and comfort level.

### ROM rating

Identify the area and surrounding muscles according to the physiotherapy protocol. Start the instrument and set it at low to medium intensity. The pressure applied is the weight of the instrument itself.

### Lymphatic drainage

Requires light, generalized physiotherapy.

Move slowly from the distal to the proximal end of the limb, i.e. in the direction of the cardiac venous return.

Reduces swelling and eliminates pain.

Restore muscle flexibility and extensibility.

Promote blood circulation and lymphatic return.

### Postoperative

The first postoperative physiotherapy treatment is based on the patient's response and the gradual increase in the duration and intensity of the treatment.

Determine the site of physiotherapy.

Ask the patient if there are any open wounds, stitches to drain the internal fixation. (If so, do not stimulate these areas)

Start the device and set it on the lowest setting.

Cover the treatment area with a towel.

Apply moderate pressure to the physiotherapy site with the instrument.

Limit physiotherapy to the area around the surgical site.

### Efficacy

Elimination of post-operative and post-injury purpura tissue

Trigger point physiotherapy.

## Safety instructions

**DANGER: To reduce the risk of electric shock, fire and personal injury, this product must be used in accordance with the following instructions**

For adult use only

- 1.Can only be used on the dry, clean skin surface of the body parts, not through the clothing gently pressed on the skin, moving, each part for about 60 seconds.
- 2.Use the device only on the soft tissues of the body without causing pain or discomfort Do not use on the head or any hard or bony part of the body.
- 3.Use only the vibration head that produces the best results  
Either controlled downward pressure or pressure applied may produce bruising. Check the treatment area frequently and stop using it as soon as there are signs of pain or discomfort Keep fingers, hair or other body parts away from the shaft and back of the vibrating head, as they may get caught Do not place any objects in the ventilation holes of the device
- 4.Do not immerse the device in water or allow water to enter the ventilation holes of the device Do not drop or improperly use the device
- 5.Check the device and electricity carefully before each use.
- 6.Do not alter or modify the device in any way.
- 7.Never leave the device running or charging unattended.

## Method of use

Local muscle relaxation

Static technique

Dynamic techniques

Pectoralis -30s

Deltoid -30s

Rotator cuff muscles -25s

Biceps -30s

Brachialis -30s

Brachioradialis -30s

Abdominals -45s

Iliopsoas muscle -25s

Quadriceps -60s

Hip abductors -60s

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Peroneus longus -30s

Peroneus brevis shortus -30s

Scapularis elevator -1min to 1.5min

Trapezius muscle -1min to 1.5min

Rhomboids -1min

Rotator cuff muscle -25s

Triceps -25s

latissimus dorsi - 45s

erector spinae -35s

Henricus medius 1min to 1.5min

Gluteus maximus 1min to 1.5min

Pear muscle 1min to 1.5min

Iliotibial bundle -25s

Biceps brachii -1min

Greater adductor muscle -25s

Semitendinosus&semimembranosus -40s

Hallux valgus muscle -20s

Flexor digitorum brevis -40s

Heel -20s

## Upper extremity operation program

Supina position

Round head vertical down

Brachioradialis 10 s – 15 s, 3 times

Radial carpal flexor 10 s -15 s, 3 times

Longissimus dorsi 12 s, 15 s

Ulnar carpal flexor 10 s -15 s, 3

Finger extensors 5 s -10 s<sup>2</sup>

Radial carpal short extensors 5 s -10 s<sup>2</sup>

Radial carpal long extensors 5 s 10 s, 2 times

Ulnar carpal extensors 5 s -10 s, 2 times

Elbow muscles 10 s - 15 s, 2 times

Biceps 15 s - 20 s, 3 times

Humerus 15 s - 20 s, 3 times

Rostral humerus 15 s - 20 s, 3 times

Triceps 15 s - 20 s, 3 times

Deltoid

-Anterior fasciculus 10 s - 20 s, 3 times

-Middle fasciculus 10 s - 20 s, 3 times

-Posterior fasciculus 10 s - 20 s, 3 times

## Lower extremity operation scheme

Prone position

Small flat head vertical down

Rouge cord muscle

Biceps femoris 20s-30sx3 times

Semilendinosus muscle 20s-30sx3 times

Semimembranosus muscle 20s-30sx3 times

Direction of movement from top to bottom

Small flal head 45° outward

Thin femoral muscle 10s-15sx2 limes

Longissimus dorsi 10s-15sx2 times

Short retraclors 10s-15sx2 times

Calf triceps

Overhead muscles 20s-30sX3 times

Flounder muscle 20s-30sX3 times

Moving direction from top lo bottom

Small round head vertical downward

Treadmill

Straightl knee position 15s-30s

Bend knee posilion 15s-30s

Supina position  
Small flat head vertical down

Quadriceps femoris  
Rectus femoris 30s-45sx2 times  
Lateral femur 30s-45sx2 times  
Medial femur 30s-45sx2 times  
Middle femur 30s-45sx2 times  
Direction of movement from top to bottom  
Small flat head vertical downward  
Suture muscle 30s-45s  
Broad fascial tensor fasciae 15s-20sx2 times  
Tibialis anterior muscle 15s-20sx2 times  
Peroneus longus 15s-20sx2 times  
Peroneus brevis short muscle 15s-20sx2 times  
Movement direction from top to bottom

Knee arthritis  
Supina position  
Small flat head down vertically  
Prone position  
Small round head vertical down

Quadriceps Femoris  
Reclus femoris 20s-45sx2 times  
lateral femur 20s-45sx2 times  
Medial femur 20s-45sx2 times  
Middle femoral muscle 20s-45sx2 times  
The direction of movement is from top to bottom  
and the patient slowly bends the knee slightly  
Small round head vertical down  
Iliopsoas edge - anterior iliacus 30s-1 min  
Reclus femoris, inguinalis 1 min  
Broad fascia tensor 15s-20s x 2 times  
Palellar ligament area 20s x 3 times  
Circle muscle 30s-45s  
Straight knee position 15s-30s  
Flexed knee position 15s-30s  
Calf triceps  
Overhead muscles 20s-30sx3s  
Halibut 20s-30sx3times

Tennis Elbow  
Supina position  
Small flat head vertical down

Anterior deltoid 10sx3 times  
Middle bundle 10sx3 times  
Posterior fasciculus 10sx3 times  
Biceps 15s-20sx3 times  
Brachialis 15s-20sx3 times  
Brachioradialis 15s-20sx3 times

Small round head vertical down  
Lateral humeral ankle tendon 30S40S  
If you are afraid of hilling the bone surface, you can squeeze the muscle with your hand and hit it downward  
Small round head downward vertically:

Flexor carpi radialis 10s-15sx3 times  
Palmaris longus 12s-15sx3 times  
Flexor carpi radialis ulnaris 10s-15sx3 times  
Finger extensors 5s-10sx2 times  
Flexor carpi radialis short 5s-10sx2 times  
Flexor carpi radialis longus 5s-10sx2 times

Neck and shoulder operation  
program Supine position  
Small flat head vertical down

Triceps 15s-20sx3reps  
Elbow muscles 10s-20sx2 reps  
Superior Irapezius bundle 15s-20sx3 reps  
Middle fascicle 15s-20sx3 reps  
Inferior fasciculus 15s-20sx3 reps  
Lalissimus dorsi 30s-1 minx2 reps

Direclion from the middle to the sidas, from top to bottom  
Small round head vertical down

Supraspinalus 10s-15sx3reps  
Scapular lift 10s-12sx2 times  
Large circular muscle 10s-12sx2 reps  
Small round muscle 10s-12sx2 times  
Infraspinalus 10s-12sx2 times  
Direclion from the middle to the sides

Prone position  
Chest, triceps

Upper arm open 90 degrees, wrist straight.  
Forearms perellel to the torso  
Slimulale lhe pectoralis mejor first 30s  
Switch to stimulating the brachialis 30s

Low back operation program  
Supina position  
Small flat head down vertically

Oblique muscle (middle and lower fibers)  
Middle fasciculus 15s-20sx3 reps  
Inferior fasciculus 15s-20sx3 reps  
Rhomboids 20s-30sx3 reps  
Vertical spinal muscles (lumbar segment) 30s-1 minx4 reps  
Iliocostalis lumborum 30s-1 minx4 reps  
Thoracolumbar fascia 1minx3 times  
Sacral joint 1minx3 times  
Small flat head vertical down  
Truncalis major 30s-1 minx6 reps  
Supine position  
Internal and external obliques 10s-15sx2 times  
Small round head vertical down erector spinae (thoracic segment) 15s-20sx2 reps  
Lumbar square 20s-30sx2 times  
Iliopsoas muscle 30s-40sx4 times

Lumbar muscle strain  
Prone position (unilateral) Small  
flat head down vertically

Gluteus maximus 30s-1 minx4 times  
Gluteus minimus and skeleton 30s-1 minx3 times  
Sacral skeletal joint 30s-1 minx2 times  
Lateral position, bend the straightened leg below by 90  
Small flat head down vertically:  
Gluteus medius and skeleton 30s-1 r  
Sacral skeletal joint 30s-1 r  
Skeletal tibial bundle 30s-1 r  
Around greater trochanter 30s-1 r  
Supine position (unilateral)  
Skeletal muscles 30s-1 minx3 times  
Slow flexion of the patient's hip  
Skeletal margin-anterior skeleton 30s-1 minx2 times  
Rectus femoris inguinalis 1 minx3 times

Sacral skeletal operation protocol  
Prone position (unilateral)  
Small flat head vertical down

Gluteus medius 30s-1 minx4 times  
Gluteus maximus 30s-1 minx4 times  
Paretic muscle 1 minx1 .5min  
Sacral joint 1 minx1 .5min

Direction: inward upward to outward downward, along the muscle

Rousseau muscle  
Biceps femoris 20s-30sx3 times  
Semitendinosus 20s-30sx3 times  
Semimembranosus 20s-30sx3 times  
From top to bottom

Plantar fasciitis in  
supina position  
Small flat head vertical down

Orbital cord muscle  
Biceps femoris 20s-30sx3 times  
Semitendinosus 220s-30sx3 times  
Semimembranosus muscle 220s-30sx3 times  
Direction of movement from top to bottom  
U-shaped head vertical downward circling muscle  
Straight knee position 215s-30s  
Flexed knee position 215s-30s  
Calf triceps  
Gastrocnemius 220s-30sx3 reps  
Hallux valgus 220s-30sx3 reps  
Four-point position with forefoot on the ground  
Small round head perpendicular to the sole of the foot  
Plantar fascia 220s-30sx3 reps  
Supina position  
Tibialis anterior 215s-20sx2 reps  
Peroneus longus 215s-20sx2 times  
Peroneus brevis 215s-20sx2 times  
Movement direction from top to bottom

**Warranty Policy**

This product complies with the warranty law in force in the country of purchase. When you have a problem that you can not solve, you must access the web [www.prixton.com](http://www.prixton.com) and click on the contact option to send us your assistance form.

These technical and electrical specifications are those of typical units and may be changed without prior notice in order to provide improved equipment.

**Simplified EU Declaration of Conformity**

We, La Trastienda Digital located in Avenida Cervantes 49, módulo 25, Basauri 48970, Bizkaia hereby declare that the following goods, imported and selling at European market:

<b>BRAND</b>	<b>MODEL</b>	<b>DESCRIPTION</b>	<b>CATEGORY</b>
PRIXTON	POWER FLEX	MASSAGE GUN	FITNESS

Are in compliance with the EC Council Directive of 2014/53/EU (EM DIRECTIVE) and harmonized standards and with the board of the 16th of April 2014.

The complete text of the UE Declaration of conformity is available anytime at this website address:

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